

DL-300 Wireless

Wireless scanning, endless possibilities



Fast & Accurate
Wireless Scanning



Outstanding Battery
Performance



Up to 8m Wireless
Range



New Software with
Advanced Features



Calibration Free



Remote Control

Fast & Accurate Wireless Scanning


DL-300 Wireless adopts cutting-edge wireless data transmission technology, allowing you to scan quickly and stably from any angle.

Outstanding Battery Performance


With up to 1 hour of continuous scanning on a single charge, you can enjoy uninterrupted and smooth scanning.

Up to 8m Wireless Range

No wires to hold you back! Enjoy the freedom of wireless movement with a range up to 8 meters.



Forma
dental supplies



New Software with Advanced Features

The brand-new software UI with built-in tutorials helps you get started quickly and easily.

Calibration Free

DL-300 Wireless is a calibration-free scanner designed to offer hassle-free scanning without the need for regular calibration adjustments. The scanner is ready whenever you need it while ensuring accurate and consistent performance.



Remote Control

Enjoy a more hygienic workflow and convenient controls at your fingertips. The DL-300 wireless allows you to effortlessly check scan data by controlling the scanner itself. Rotate, zoom in, rescan, or switch between interfaces with a simple gesture without interrupting your flow or concentration.

What's in the box...

Intraoral Camera



* The scanner & battery are packed separately

Power Adapter



Scanner Tips



* 1 Tip installed on the scanner already

User Manual



Camera Holder



Batteries
(3 included)



Dongle



LAUNCA


Forma
dental supplies

