

EzScan™



Expand Your Digital Workflow



Simple and Intuitive Interface for an Improved Scanning Experience

- Clear and intuitive interface for ease of use
- Essential functionality for many treatment plans



Case Sharing

- EzScan Cloud is an easy-to-use file-sharing module for sharing digital impressions with the lab of your choice.



Quick Scan Time for Workflow Efficiency

- Full mouth scans in under 2 minutes from start to finish




Open Platform with No Continued Fees

- Easy to export to your treatment planning software or to send it to the lab of your choice
- Maximize your ROI with no continued fees



Bring simplicity to your workday with an ergonomically designed lightweight scanner

The EzScan is a digital 3D imaging solution designed to bring simplicity to your work. It is presented as a turn-key solution that takes under 10 minutes to unbox and setup. The solution includes a pre-calibrated scanner, three autoclavable tips, and an optimized acquisition PC with the EzScan software pre-installed to ensure that you spend less time installing and more time scanning. The EzScan-i features an intuitive, user-friendly interface that makes it easier than ever to meet all of your restorative, orthodontic, and implant needs.



360°
rotation

59 sec.
for a full arch scan

up to
250 times
autoclavable tips

150g
small & lightweight

SPECIFICATION

FOV	12 x 14 mm	Features	True Color / Powder Free
Scan Time (Full Arch)	Less than 2 min.	File Output	STL, PLY, OBJ, DWOS
STL Converting Time (Full Arch)	Less than 1 min.	System	Open
Weight	150 g	Connector	USB Type-A