

Medit i600

Smart versatile scanning

Meet the Medit i600, where style meets versatility. Leading in dental technology, it offers smart scanning with ease of use for any procedure, making it essential for dental practices.



Inside the Precision: Medit i600

Discover Medit i600's core performance and features.



35 FPS
Stable scanning

Reversible Tip
180° flexibility

High Precision
10.9 μm ± 0.98

Scan Button
Easy start-up

245g
Lightweight

Plug and Scan
Easy connection

Lighting the way in dental scanning

The Medit i600 shines a light on new frontiers in dental care, standing as a beacon of innovation. It paves the way for a more efficient, precise, and insightful approach to dental scanning, redefining the standards of modern dentistry.





Forma
dental supplies

Durably lightweight, 245g

Built to last, easy to handle

Medit i600 is a sleek blue scanner weighing 245 grams, both lightweight and powerful.



Agile and blue

Reduce hand strains

Enhance your comfort and reduce hand strain during dental procedures.



Zero hassle setup

Just plug and scan

Connect directly with USB Type-C, no power hubs needed.



Scan swiftly with confidence

Faster scans with both accuracy and assurance

Deliver precise, detailed outcomes confidently with reliable data. Enjoy quick scanning without compromising your patients' trusted quality.



Up to 35 fps

Smooth scanning experience

Achieve quick, precise scanning with a steady 35 fps, matching the [Medit i700](#).



High accuracy

Capture precise scan data

Achieve clear, detailed scans with 10.9 μm accuracy, ensuring top optical quality.



180° reversible tips

Improved convenience

Easily adapt to various oral conditions with the reversible tip.



Maximum versatility

Smart scanning for every need


Handle full-arch scans to intricate restorative work with smart performance.



Better hardware with software updates

Innovate with Medit, where updates mean upgrades

Software updates transform your scanners to perform better.



Transform your dental practice today

Request a demo and experience it firsthand

See our intraoral scanner in action. Request a demo to transform your dental practice today.

