

PLANMECA

Planmeca Pro50™ Chair

The patient chair for premium comfort



Planmeca Pro50™ Chair is a great combination of functionality, flexibility and future-proof details, all wrapped up in a beautiful, timeless design. It is an ideal option for those looking for a sturdy yet elegant patient chair – offering ultimate comfort for both the dental team and patients. Have a seat!

PLANMECA

Silent and sturdy

Planmeca Pro50™ Chair is a modern and timeless patient chair featuring all the lightweight and beautiful design details that Planmeca products are known for. The medical grade motors inside the chair make its movements smooth and silent, ensuring a quiet working environment for the dental team and a relaxing treatment experience for the patient. The chair is built from durable aluminium and its upholstery comes in a wide range of fun and appealing colours.

First-class ergonomics

Planmeca Pro50 Chair offers a wide height range – the patient seat height can be varied between 390 mm and 890 mm. This allows dental professionals of all sizes to work sitting or standing. The small chair base leaves plenty of room for two pairs of legs as well as the foot control and other devices. The slim and narrow backrest allows the dental team to get close to the patient – always in the most ergonomic working position.



PLANMECA

Premium patient comfort

Planmecca Pro50 Chair has a folding leg rest that provides easy entry and exit and facilitates face-to-face positioning – a feature that is greatly appreciated in doctor-patient consultations. As patients come in different shapes and sizes, the maximum patient weight for the chair is 185 kg. The Ultra Relax™ upholstery adapts perfectly to the patient's body shape – allowing them to stay relaxed even during the longest treatment sessions.



A smooth and flexible workflow

For those working with a 12 o'clock cabinet or tray, the chair's lifting mechanism ensures that the distance from the headrest to the cabinet remains the same at every height and the operating light doesn't need to be readjusted as the chair height changes. The tilted design of the chair ensures optimal patient positioning when moving from a sitting to a supine position. What's more, the chair has an optional mobile base, making it easy to move the chair to a new treatment room when necessary.

Timeless and future-proof


Just like all Planmecca products, Planmecca Pro50 Chair is guaranteed to have a long life span. This, combined with its timeless design and sturdy metal construction, makes it a truly future-proof investment.




Other Planmeca dental unit features and specialities



Build your own unit



Dental unit accessories



Planmeca Emerald® S