

# PLANMECA

## Planmeca ProScanner® 2.0

User-friendly and reliable imaging plate scanner



The Planmeca ProScanner® 2.0 imaging plate scanner is a high quality option for fast and dependable intraoral imaging. With intelligent design details and outstanding durability, the scanner supports everyday tasks at a dental clinic – providing reliability that is hard to match.

## Reliable and dependable

Planmeca ProScanner 2.0 has been designed to be a seamless part of a clinic's workflow. It is a true workhorse imaging plate scanner that is maintenance-free and extremely durable. The compact scanner fits into any room, helping all clinics prevent downtime and ensure maximum efficiency.



## Smart imaging plates

Planmeca ProScanner 2.0 utilises RFID technology for fast scanning and identification. The imaging plate scanner's high-quality plates can be used numerous times and the integrated eraser instantly prepares them for further imaging. The flexible imaging plates are very comfortable for patients and available in the most common sizes – 0, 1, and 2.



## Easy to use

Planmeca ProScanner 2.0 significantly speeds up the intraoral imaging workflow with its fast scanning process. Just insert the imaging plate, scan it and diagnose the captured image. If the wrong side of the plate is accidentally exposed, it is easy to notice immediately from the captured image.



## Planmeca Romexis® software integration

Each imaging plate is equipped with a digital serial number that allows sorting, viewing and comparing images for quality control, as well as counting exposures by using the serial number of each plate in the Planmeca Romexis® software.

