

# PLANMECA



## Planmeca ProSensor® HD

Intraoral sensor with three different sizes

The Planmeca ProSensor® HD intraoral sensor offers a unique combination of unparalleled image quality, high-end patient-centred design and advanced usability. The digital sensor has been designed to last and provides durability that can be counted on.

### Cutting-edge image quality

Supporting detailed diagnostics, Planmeca ProSensor HD captures sharp, low-noise and high-contrast images. With a true resolution of over 20 lp/mm, the intraoral sensor delivers HD image quality and its wide dynamic range ensures successful results each and every time.



## Patient-centred design

To fulfill all intraoral imaging needs, Planmeca ProSensor HD's sensors come in three different sizes. Their rounded edges make procedures comfortable for patients, with outstanding images ready in a matter of seconds.



## High-end usability

Planmeca ProSensor HD's magnetic connector makes the intraoral sensor simple to attach, with its white colour enhancing visibility to further ease positioning. The sensor's elegant control box is equipped with a colour-coded LED light that provides instant visual feedback of the imaging procedure, while hermetically sealed housing ensures effective infection control.

## Designed to last

Intraoral sensors are naturally subject to wear and tear in daily use, but Planmeca ProSensor HD has been built to last. The sensor's carbon fibre plate prevents bite marks, its shock protection layer guards against minor damage and its strengthened cable further improves durability.



## Intraoral X-ray unit integration

The seamless integration of the Planmeca ProSensor HD intraoral sensor and the Planmeca ProX™ intraoral X-ray unit ensures efficiency at every step. Due to its smooth magnetic connector, the sensor can be easily attached using just one hand.

