


PLANMECA

Planmeca Solanna®

Flexible dental unit operating light



The Planmeca Solanna® dental unit operating light provides perfect visibility over the entire treatment area – with the brightness and colour temperature of your choice. The light's beautiful and effortless movements allow smooth positioning in all situations.

Uniform light with less shadows

Planmeca Solanna's advanced LED technology and unique mirror optics guarantee uniform light intensity and minimised shadows, regardless of the distance from the treatment area.



PLANMECA

Clear visibility to all corners

The dental unit operating light's large illumination field provides excellent visibility to every corner of the treatment area, minimising the need to adjust the light head during treatments. The beam's sharp edges make the light eye-friendly for patients.

Easy positioning

Planmeca Solanna's smoothly moving light head caters to the needs of users of all sizes – sitting or standing. The light head's three-axis movement makes it easy and effortless to position and the exceptionally long arm enables optimal working conditions for operations of all kind.



Adjustable intensity and colour

Planmeca Solanna's light intensity can be adjusted flexibly between 8,000 and 50,000 lux, with three colour temperature options – warm, neutral and cool. The light's composite mode enables safe use of composite materials and prevents pre-polymerisation.

PLANMECA



Different mounting possibilities

Planmeca Solanna is a flexible dental unit operating light that offers multiple mounting options. The light can be attached to either a dental unit, the ceiling or a wall.



Effective infection control

You can wash Planmeca Solanna's detachable and autoclavable handles in a washer-disinfector. The light's seamless one-piece aluminium back cover and dustproof light capsule further guarantee optimal infection control.

PLANMECA

Usability at its finest

Planmeca Solanna can be controlled in several ways – from its innovative handles, the dental unit's user interface or foot control, or with defined hand movements. Connecting the light to the Planmeca Romexis® software provides reliable usage information.



Other Planmeca dental unit features and specialities



Build your own unit



Dental unit accessories



Planmeca Emerald® S